Mediblink[®]
medical devices

EN Mediblink Blood Pressure Monitor M520 AFib

INSTRUCTIONS FOR USE

Please read the instructions for use carefully before using the product

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Pažljivo pročitajte upute prije upotrebe uređaja

TABLE OF CONTENTS

| | |
|---|-----------|
| 1. Introduction | 6 |
| 1.1 Features | 6 |
| 1.2 Important information about self-measurement | 6 |
| 2. Important Information of Blood-Pressure Measurement | 7 |
| 2.1 How does high / low blood-pressure arise? | 7 |
| 2.2 Which values are normal? | 7 |
| 3. The Various Components of the Blood-Pressure Monitor | 9 |
| 4. Putting the Blood-Pressure Monitor Into Operation | 10 |
| 4.1 Inserting the batteries | 10 |
| 4.2 User selection and setting the time / date | 10 |
| 5. Carrying Out a Measurement | 11 |
| 5.1 Before the measurement | 11 |
| 5.2 Common sources of error | 11 |
| 5.3 Fitting the cuff | 12 |
| 5.4 Measuring procedure | 13 |
| 5.5 Discontinuing a measurement | 14 |
| 5.6 Memory – storage and recall of the measurements | 15 |
| 5.7 Memory – cancellation of all measurements | 15 |
| 6. Appearance of the Pulse Arrhythmia Indicator for Early Detection | 16 |
| 7. Appearance of the Atrial Fibrillation Indicator for Early Detection | 16 |
| 8. Error Messages / Malfunctions | 17 |
| 9. Care and Maintenance, Recalibration | 18 |
| 10. Battery Life | 19 |
| 11. Safety, Care and Disposal | 19 |
| 12. Reference to Standards | 20 |
| 13. Remark | 20 |
| 14. Technical Specifications | 21 |
| 15. Manufacturer's Declaration | 22 |
| 16. Warranty | 24 |

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1. Introduction

1.1 Features

The Mediblink Blood-pressure monitor M520 AFib with technology AFib (Atrial Fibrillation), MDI (Measuring During Inflation), integrated time and date display, is a fully automatic, digital blood-pressure measuring device for use on the arm, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measurement.

The device offers very high and clinical tested measurement accuracy and has been designed to provide a maximum of user-friendliness.

The device is intended for self-use in home.

Before using, please read through this instruction manual carefully and then keep it in a safe place.

For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

Warning!

1.2 Important information about self-measurement

- Substitution of a different component might result in measurement error.
- Cuff is replaceable only by an original.
- Do not use with neonatal patients.
- It will cause harmful injury to the patient or affect the blood pressure due to connection tubing kinking.
- Too frequent measurements can cause injury to the patient due to blood flow interference.
- The application of the cuff over a wound can cause further injury.
- The application of the cuff and its pressurization on any limb where intravascular access or therapy, or an arteriovenous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury to the patient.
- Do not let the cuff and its pressurization on the arm on the side of a mastectomy.
- The need to check that operation of the automated sphygmomanometer does not result in prolonged impairment of patient blood circulation.
- Not intended to be used together with HF surgical equipment.
- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.

- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. Important Information of Blood-Pressure Measurement

2.1 How does high / low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory center, and adapted to the respective situation by way of feedback via the nervous system. To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls. The level of arterial blood-pressure changes periodically during the heart activity: During the "blood ejection" (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's "rest period" (Diastole) minimal (diastolic blood-pressure value). The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2 Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and / or the systolic blood-pressure is over 160 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

Should the systolic blood-pressure values lie between 140 mmHg and 160 mmHg and / or the diastolic blood-pressure values lie between 90 mmHg and 100 mmHg, likewise, please consult your doctor. Furthermore, regular self-checks will be necessary.

With blood-pressure values that are too low, i.e. systolic values under 100 mmHg and / or diastolic values under 60 mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately. If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.**


Table for classifying blood-pressure values (unit: mmHg) according to World Health Organization:

| Range & broadcasting | Systolic Blood-pressure | Diastolic Blood-pressure | Measures |
|----------------------|-------------------------|--------------------------|-------------------------------|
| Optimal | Between 100 and 120 | Between 60 and 80 | Self-check |
| Normal | Between 120 and 129 | Between 80 and 84 | Self-check |
| High to normal | Between 130 and 139 | Between 85 and 89 | Consult your doctor |
| Slight hypertension | Between 140 and 159 | Between 90 and 99 | Seek medical advice |
| Medium hypertension | Between 160 and 179 | Between 100 and 109 | Seek medical advice |
| Strong hypertension | Higher than 180 | Higher than 110 | Urgently seek medical advice! |

Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120mmHg require immediate medical treatment.

3. The Various Components of the Blood-Pressure Monitor

4. Putting the Blood-Pressure Monitor Into Operation

4.1 Inserting the batteries

- Insert the batteries (4 x size AA 1.5V), thereby observing the indicated polarity.
- If the battery warning “” icon appears in the display, the batteries remain 20% power to warn user the batteries will be run out.
- If the battery warning “” icon appears in the display, the batteries are empty and must be replaced by new ones


Warning!

- After the battery warning “” icon appears, the device is blocked until the batteries have been replaced.
- Please use “AA” Long-Life or Alkaline 1.5V Batteries.
- The use of 1.2V Accumulators is not recommended.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

4.2 User selection and setting the time / date

User selection: This advanced blood pressure monitor allows you to track blood pressure readings for 2 individuals independently.

- Before measurement, make sure you set the unit for the intended user. The unit can track results for 2 individuals (User 1, User 2).
- Press the Memory button “” for at least 3 seconds. The display now indicates the set user, during which the set user blink, to confirm, press ON / OFF button.
- Press the MEMORY button to select User.
- We suggest the first person to take their pressure to be User 1.

Setting the time, date

This blood-pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood-pressure values are stored, but also the exact moment of the measurement.

How to enter to the setting mode:

After batteries are inserted, the device will enter the setting mode automatically.

The setting procedure as follows:

It is essential to set date and time. Otherwise, you will not be able to save your measured value correctly with a date and time.

To use “**M**” button to be adjustment button.

If you press and hold “**M**” button, you can set the values more quickly.

1. Year setting: the year flashes on the display: choose the desired year with “**M**” button and confirm with “**⏻**” (ON / OFF) button.
2. Month setting: the month flashes on the display: choose the desired month with “**M**” button and confirm with “**⏻**” (ON / OFF) button.
3. Day setting: the day flashes on the display: choose the desired day with “**M**” button and confirm with “**⏻**” (ON / OFF) button.
4. Hour setting: the hour flashes on the display: choose the desired hour with “**M**” button and confirm with “**⏻**” (ON / OFF) button.
5. Minute setting: the minute flashes on the display: choose the desired unit with “**M**” button and confirm with “**⏻**” (ON / OFF) button.

Please be noticed, the clock will begin to run from 2016-01-01 01:01, if no key is pressed within 20 seconds.

5. Carrying Out a Measurement

5.1 Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about ten minutes before the measurement.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

5.2 Common sources of error

Note: Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- The performance of the automated sphygmomanometer can be affected by extremes of temperature, humidity and altitude.
- Avoid compression or restriction of the connection tubing.
- A loose cuff causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results. Correctly executed blood-pressure measurements should therefore first be repeated after a 5 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away (after at least 3 minutes).

5.3 Fitting the cuff

Insert air connector into air outlet shown in left photo and please make sure the fitting of the air connector completely and properly to avoid air leakage.

- The distance between the edge of cuff and the elbow should be approx. 1-2cm.
- Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight, whereby no space should remain between the cuff and the arm.

- Lay the arm on a table, with the palm upwards. Support the arm a little with a rest (cushion), so that the cuff rests at about the same height as the heart. Take care, that the cuff lies free. Remain so for 2 minutes sitting quietly, before beginning with the measurement.


- Let legs uncrossed, feet flat on the floor, back and arm supported.

5.4 Measuring procedure

5.4.1 Measuring in standard mode

In this mode the device has arrhythmia detection but no AFib detection.

After the cuff has been appropriately positioned, the measurement can begin:

- 1) Press the **ON/AFib** button, the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- 2) **Cuff fitting detection:** the icon  will appear and blink during measuring, if cuff was fit too loose. The icon  will appear during measuring, if cuff was fit well.
- 3) **Arm movement detection during measuring:** the icon  will appear, if a movement was detected which may influence accuracy. Due to the movement not too serious, the measuring can be continuous (if the movement is too serious, **Err2** displayed).
- 4) After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink
- 5) When the measurement has been concluded, the measured systolic and diastolic blood-pressure values as well as the pulse frequency are now displayed.
Example (Fig.): Systole 126, Diastole 85, Pulse 78
- 6) The measurement results are displayed, until you switch the device off or the device switches automatically off, to save the batteries.

5.4.2 Measuring in AFib Mode (2-measurement mode)

In AFib mode, 2 measurements are automatically taken in succession and the result is then automatically analyzed and displayed. Because blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement.

- Hold the **ON/AFib** button for more than two seconds until the symbol  appears on the display.
- The middle, left hand section of the display shows a 1, 2 to indicate which of the 2 measurements is currently being taken.
- There is a break of 15 seconds between the measurements.
A countdown indicates the remaining time.
- The individual results are not displayed. Your blood pressure will only be displayed after all 2 measurements are taken.
- Do not remove the cuff between measurements.

- If one of the individual measurements was questionable, a third one is automatically taken.

In the measuring:

After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink.


Measured result:

The measured systolic and diastolic blood-pressure values as well as the pulse are now displayed.

Example 1:

Systole 120, Diastole 80, Pulse 80, AFib detected.

Icon of Arrhythmia  and AFib **AF** detected;
Arm movement detected and cuff fit too loose detected

Example 2:

Systole 120, Diastole 80, Pulse 80, arrhythmia  detected, but no AFib detected.

Arm movement detected and cuff fit well.

The measurement results are displayed until you switch the device off. If no button is pressed for 3 minutes, the device switches automatically off.

5.5 Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the "ON / OFF" power button can be pressed at any time. The device then immediately lowers the cuff-pressure automatically.

5.6 Memory – storage and recall of the measurements


The blood-pressure monitor automatically stores each of the last 120 measurement values. By pressing the MEMORY button, an average value of the last 3 measurements as well as the last measurement and the further last 120 measurements (M 119, M 118 ... M 1) can be displayed one after the other.

M A
Average value
of the last 3
measurements

M 3
Measurements
value of
memory 3


M 2
Measurements
value of
memory 2

M 1
Measurements
value of
memory 1

5.7 Memory– cancellation of all measurements

Before you delete all readings stored in the memory, make sure you will not need refer to the readings at a later date. Keeping a written record is prudent and may provide additional information for your doctor's visit. In order to delete all stored readings, depress the MEMORY button for at least 5 seconds, the display will show the symbol "CL" and then release the button. To permanently clear the memory, press the MEMORY button while "CL" is flashing. To indicate deletion of stored readings.

6. Appearance of the Pulse Arrhythmia Indicator for Early Detection

This symbol “” indicates that certain pulse irregularities were detected during the measurement.

In this case, the result may deviate from your normal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor.

Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Arrhythmia indicator

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested. The arrhythmia symbol is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

7. Appearance of the Atrial Fibrillation Indicator for Early Detection

This device is able to detect atrial fibrillation (AFIB). This icon “” indicates that atrial fibrillation was detected during the measurement. If the AFib symbol appears after having performed a full blood pressure measurement episode (double measurements), you are advised to wait for one hour and perform another measurement episode (double measurements). If the AFib symbol appears again, then you are advised to visit your doctor. If after repeated measurement the AFib symbol is no longer displayed there is no cause for concern. In such case it is recommended to measure again the next day. Keep the arm still during measuring to avoid false readings. This device may not detect atrial fibrillation in people with pacemakers or defibrillators.

8. Error Messages / Malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.

| Error No. | Possible cause(s) |
|-----------|---|
| ERR 1 | No pulse has been detected. |
| ERR 2 | Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the Measurement (Artefact). |
| ERR 3 | The inflation of the cuff takes too long. The cuff is not correctly seated. |
| ERR 5 | The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take other reading following directions carefully. Contact you doctor if you continue to get unusual readings. |
| ERR 8 | The pressure in the cuff is too high (over 290 mmHg). |

Further Information – The level of blood-pressure is subject to fluctuations even with healthy people. Important thereby is, that comparable measurements always require the same conditions (Quiet conditions)! If, in spite of observing all these factors, the fluctuations are larger than 15 mmHg, and / or you hear irregular pulse tones on several occasions, please consult your doctor. For licensing, the device has been subjected to strict clinical tests, by which the computer program used to measure the blood-pressure values was tested by experienced specialist doctors in Germany. The same computer program is used in every individual device, and has thus also been clinically tested. The manufacture of the devices takes place according to the terms of the European standard for blood-pressure measuring devices (see technical data) you must consult your specialist dealer or chemist if there are technical problems with the blood-pressure instrument. Never attempt to repair the instrument yourself! Any unauthorized opening of the instrument invalidates all warrantee claims!

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

| Malfunction | Remedy |
|---|---|
| The display remains empty when the instrument is switched on although the batteries are in place. | <ol style="list-style-type: none"> 1. Check batteries for correct polarity and if necessary insert correctly. 2. If the display is unusual, re-insert batteries or exchange them. |
| The device frequently fails to measure the blood pressure values, or the values measured are too low (too high). | <ol style="list-style-type: none"> 1. Check the positioning of the cuff. 2. Measure the blood-pressure again in peace and quiet under observance of the details made under point 5. |
| Every measurement produces a different value although the instrument functions normally and the values displayed are normal | <ol style="list-style-type: none"> 1. Please read the following information and the points listed under "Common" sources of error. Repeat the measurement. Please note: Blood pressure fluctuates continually so successive measurements will show some variability. |
| Blood pressure measured differs from those values measured by the doctor. | <ol style="list-style-type: none"> 1. Record the daily development of the values and consult your doctor. Please note: Individuals visiting their doctor frequently experience anxiety which can result in a higher reading at the doctor than obtained at home under resting conditions. |

9. Care and Maintenance, Recalibration

- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- Never open the device! Otherwise the manufacturer calibration becomes invalid!

10. Battery Life

1000 times measurement with 4- size "AA" alkaline Batteries.

11. Safety, Care and Disposal

Safety and protection

- This instrument may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for the damage caused by incorrect application.
- This instrument comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the "Technical specifications" section!
- Protect it from water and moisture, extreme temperatures, impact and dropping, contamination and dust, direct sunlight, heat and cold.
- The cuffs are sensitive and must be handled with care.
- Only pump up the cuff once fitted.
- Do not use the instrument close to strong electromagnetic fields such as mobile telephones or radio installations.
- Do not use the instrument if you think it is damaged or notice anything unusual.
- If the instrument is not going to be used for a prolonged period the batteries should be removed.
-  Read the additional safety instructions in the individual sections of this booklet. Ensure that children do not use the instrument unsupervised: some parts are small enough to be swallowed.
- Must use the recognized accessories, detachable parts and materials, if the use of other parts or materials can degrade minimum safety.
- A warning to remove primary batteries if the instruments is not likely to be used for sometime.

Instrument care

Clean the instrument only with a soft, dry cloth.

Disposal

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

12. Reference to Standards

Device standards: Device corresponds to the requirements of the European standard for non-invasive blood- pressure monitor.

Standards

EN IEC 80601-2-30, EN ISO 81060-2, EN 60601-1, EN 60601-1-6, EN 62366-1, EN 62304, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-23, ISO 14971, EN 60601-1-2, EN 60601-1-11

The Device meets the requirements of Directive 93/42/EEC for Medical Devices Class IIa.

13. Remark

| | | | |
|---|--|---|---|
|  | Electrical devices are recyclable material and should not be disposed of with household waste after use! |  | Serial number |
|  | Manufacturer |  | Medical Device |
|  | Inapplicable for baby |  | Unique Device Identifier |
|  | Cuff Connector |  | Number of products in one packaging |
|  | Warning! Consult accompanying documents |  | Before use, read the instructions. Electronic instructions for use: http://www.mediblink.com/ff/m520.pdf |
|  | Product reference number |  | Type BF equipment |
|  | Model number |  | Degrees of protection |
|  | Batch number* |  | Authorised Representative in EU |

*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

14. Technical Specifications

| | |
|---|--|
| Model/REF: | BP168A/M520 |
| Measurement Procedure: | Oscillometric, corresponding to Korotkoff method: |
| Phase I: | Systolic, Phase V: diastolic |
| Display: | Digital display |
| Measuring range: | Pressure: SYS: 60 to 255 mmHg; DIA: 30 to 199 mmHg (in 1 mmHg increment) Pulse: 40 to 199 beat / minute |
| Static accuracy: | Pressure: ± 3 mmHg Pulse: ± 5 % of reading |
| Measuring resolution: | 1 mmHg |
| Inflation: | Automatic inflation by internal pump |
| Memory function: | 2 x 120 memories for 2 users (SYS, DIA, Pulse) |
| Decompression: | Constant exhaust valve system |
| Power source 1 (Battery): | 4- size "AA" alkaline Batteries |
| Power source 2 (USB-C Adapter): | Input : 100-240V ~ 50 / 60 Hz 0.2A Output : 5V / 1A, 0.5W |
| Operation temperature: | 5 ~ 40°C / 41 ~ 104°F |
| Operation humidity: | 15 % ~ 80 % RH maximum |
| Storage temperature: | -20 ~ +55°C / -4 ~ +131°F |
| Storage humidity: | 10 % ~ 95 % RH maximum |
| Dimensions: | 120 × 80 × 38 ±1.0 mm |
| Weight: | 510 g ± 5g (including batteries and cuff) |
| Cuff pressure display range: | 0 ~ 290 mmHg |
| Electrical shock protection: | Internal power unit |
| Safety classifications: | Type BF equipment  |
| Mode of operation: | Continuous operation |
| Protection against ingress of water: | IP20 |
| Accessories: | M-L size Cuff, 4 "AA" batteries, storage bag, instruction manual, USB-C Adapter |

15. Manufacturer's Declaration

M520 AFib with MDI is intended for use in the electromagnetic environment specified below. The customer or the user of M520 AFib with MDI should assure that it is used in such an environment.

Electromagnetic Emissions: (IEC60601-1-2)

| Emission Test | Compliance | Electromagnetic Environment |
|--|----------------|--|
| RF emission CISPR 11 | Group 1 | M520 with MDI uses RF energy only for internal functions. Therefore, this RF emission is extremely weak and there is little chance of it creating any kind of interference whatsoever with nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | M520 with MDI is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations / flicker IEC 61000-3-3 | Not applicable | |

Electromagnetic Immunity: (IEC60601-1-2)

| Immunity test | IEC60601-1-2 test level | Compliance level | Electromagnetic environment – guidance |
|---|--|----------------------------|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 % |
| Electric fast transient / burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input / output lines | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ±1 kV differential mode ±2 kV common mode | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5 % U_T (95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (95 % dip in U_T) for 5 sec. | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. If the user of the upper arm style requires continued operation during power mains interruptions, it is recommended that M520 with MDI be powered from an uninterruptible power supply or a battery. |
| Power frequency (50 / 60 Hz) magnetic field IEC 61000-4-8 | 3 A / m | Not applicable | Not applicable |

Note: U_T is the a.c. mains voltage prior to application of the test level.

| Immunity test | IEC60601-1-2 test level | IEC60601-1-2 test level | Electromagnetic environment – guidance |
|----------------------------|--|-------------------------|--|
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz 80% AM (2Hz) | 3 Vrms | Portable and mobile RF communications equipment should be used no closer to any part of M520 with MDI, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommend separation distance 3 V $d = 1.2 \times P^{1/2}$ 80 MHz to 800 MHz $d = 2.3 \times P^{1/2}$ MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ³ , should be less than the compliance level in each frequency range ³ . Interference may occur in the vicinity of equipment marked with the following symbol:  |
| Radiated RF IEC 61000-4-3 | 3 Vrms 80 MHz to 2.5 GHz 80% AM (2Hz) | 3 V / m | |

Note1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

*Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the M520 with MDI is used exceeds the applicable RF compliance level above, M520 with MDI should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating M520 with MDI.

*Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V / m.

Recommended Separation Distances:

Recommended separation distance between portable and mobile RF communications equipment and M520 with MDI.

M520 with MDI is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of M520 with MDI can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and M520 with MDI as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter m | | |
|---|---|---|--|
| | 150 kHz to 80 MHz $d = 1.2 \times P^{1/2}$ | 80 to 800 MHz $d = 1.2 \times P^{1/2}$ | 800 MHz to 2.5 GHz $d = 2.3 \times P^{1/2}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies

Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Manufacturer Information:

Shenzhen Combei Technology Co.,Ltd.
11-5B, No.105, Huanguan South Road,
Dahe Community, Guanhu Street, Longhua District,
Shenzhen, 518110 Guangdong, P.R. China
kevin.fong@combei.cn

Manufactured for
(Importer for EU & Distributor):
Mediblink d.o.o.
Gubčeva cesta 19
8210 Trebnje
Slovenia
info@mediblink.com
www.mediblink.com

MedNet EC-REP GmbH

Borkstrasse 10
48163 Münster
Germany

Instructions for use, version No:
M520-2019-V03
Issue date:
20.12.2024
Date of last change:
14.9.2025

4.0, 14.9.2025

0197

16. Warranty

Product: Mediblink Pressure Monitor M520 Afib

Manufactured for (Importer for EU & Distributor): Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

Sellers name, address, signature and stamp*:

Date of extradition / sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

WARRANTY TERMS

Dear customers!

The warranty period is **10 years** and starts from the day of product purchase. The warranty on cuff is **1 year**. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product can not be repaired or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product can not be repaired in 45 days, the product will be replaced with a new one. In case that the product can not be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and / or patient is established.

Replacement parts can be ordered at www.mediblink.com.

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